



**DIN 02247354** 

### **Product Description**

DIN-Registered Super Sani-Cloth® Wipe is a premoistened nonwoven durable wipe containing a quaternary ammonium chloride/ alcohol based solution. Recommended for use in hospitals and other critical care areas where the reduction of cross-contamination between treated surfaces is required. Some organisms are removed from the surface by thoroughly wiping the surface with the wipe. Most remaining organisms are killed within two (2) minutes by exposure to the liquid in the wipe.

### **Chemical Composition**

### **Active Ingredients:**

n-Alkyl (68% C <sub>12</sub> , 32% C <sub>14</sub> ) dimethyl ethylbenzyl ammonium chlorides	0.25%
n-Alkyl (60% $C_{14^{\prime}}$ 30% $C_{16^{\prime}}$ 5% $C_{12^{\prime}}$ 5% $C_{18}$ ) dimethyl benzyl ammonium chlorides	0.25%
Isopropyl Alcohol	55.00%
Other ingredients	44.50%
TOTAL	100.00%
(Does not include the weight of the wipe)	





### Efficacy

Bacterial Organism Efficacy Multi-Drug Resistant Bacteria:

Acinetobacter baumannii [ATCC 19606]

Enterobacter cloacae NDM-1 positive [CDC 1000654] ESBL Producing Escherichia coli (E. coli) [ATCC BAA-196]

Methicillin Resistant *Staphylococcus aureus* (MRSA) [ATCC 33592] Vancomycin Resistant *Enterococcus faecalis* (VRE) [ATCC 51299]

Test Method Used: Modified AOAC Germicidal Spray Method for Hard Surface Disinfection

Organic Soil Load: 5% Horse Serum or 5% Fetal Bovine Serum

Exposure Time: 2 minutes at 68-77°F

Incubation: 48 hours +/- 2 hours at 95-98.6°F

Results: No growth observed

Bacteria: Bordetella pertussis [ATCC 12743]

Burkholderia cepacia [ATCC 25416] Campylobacter jejuni [ATCC 29428] Escherichia coli (E. coli) [ATCC 11229] Escherichia coli O157:H7 [ATCC 35150] Klebsiella pneumonia [ATCC 4352] Pseudomonas aeruginosa [ATCC 15442] Salmonella enterica [ATCC 10708] Staphylococcus aureus [ATCC 6538]

Test Method Used: Modified AOAC Germicidal Spray Method for Hard Surface Disinfection

Organic Soil Load: 5% Horse Serum or 5% Fetal Bovine Serum

Exposure Time: 2 minutes at 66-77°F

Incubation: 48 hours +/- 2 hours to 6 days at 86-98.6°F

Results: No growth observed

Mycobacterium Bovis - BCG (TB): Mycobacterium bovis BCG (Tuberculosis) [ATCC 35743]

Test Method Used: Quantitative Tuberculocidal Suspension Test

Organic Soil Load: 5% Horse Serum
Exposure Time: 1 minute at 68°F
Incubation: 21 days at 98.6°F
Results: No growth observed

**Viral Organism Efficacy** 

Organic soil load:

Enveloped Viruses: Herpes Simplex type 2 [ATCC VR-734]

Human Coronavirus [ATCC VR-740] Strain 229E Influenza A2/Hong Kong (Hong Kong flu) Influenza A (H1N1) virus [ATCC VR-98]

Measles virus [Strain: Edmonston] [ATCC VR-24]

Vaccinia virus [ATCC VR-1354]

Test Method Used: This test was conducted according to U.S. Environmental Protection Agency guidelines in effect at the

time of test for determining the virucidal efficacy of disinfectants intended for use on dry inanimate surfaces.

Organic soil load: 5% fetal bovine serum. Exposure Time: 2 minutes at 68°F

Results: Virucidal according to the criteria established by the U.S. Environmental Protection Agency for

registration and labeling of a disinfectant product as a virucide.

Respiratory Syncytial virus (RSV)

Test Method Used: This test was conducted according to U.S. Environmental Protection Agency guidelines in effect at the

time of test for determining virucidal efficacy of disinfectants intended for use on dry inanimate surfaces.

5% fetal bovine serum.

Exposure Time: 2 minutes at room temperature (68°-77°F)

Results: Virucidal against Respiratory Syncytial virus (RSV) according to the criteria established by the U.S.

Environmental Protection Agency for registration and labeling of a disinfectant product as a virucide.



### Efficacy

**Viral Organism Efficacy** 

Non-enveloped viruses: Adenovirus type 5 [ATCC VR-5]

Feline Calicivirus (Surrogate for Human Norovirus) [ATCC VR-782]

Rhinovirus [ATCC VR-1110] Rotavirus Strain WA

Test Method Used: This test was conducted according to U.S. Environmental Protection Agency guidelines in effect at the time of

test for determining the virucidal efficacy of disinfectants intended for use on dry inanimate surfaces.

Organic soil load: 5% fetal bovine serum. Exposure Time: 2 minutes at 68°F

Results: Virucidal according to the criteria established by the U.S. Environmental Protection Agency for

registration and labeling of a disinfectant product as a virucide.

Bloodborne Pathogens: Hepatitis B virus (HBV) - Duck HBV

Hepatitis C virus Human (HCV) - Bovine Diarrhea virus

Test Method Used: This test was conducted according to U.S. Environmental Protection Agency guidelines in effect at the time of

test for determining virucidal efficacy of disinfectants intended for use on dry inanimate surfaces.

Organic Soil Load: Hepatitis B virus (HBV) 100% duck serum.

Hepatitis C virus (HCV) 5% horse serum

Exposure Time: 2 minutes at room temperature (68°-77°F)

Results: Virucidal against Hepatitis B and Hepatitis C virus according to the criteria established by the U.S.

Environmental Protection Agency for registration and labeling of a disinfectant product as a virucide.

HIV-1 (AIDS VIRUS)

Test Method Used: This test was conducted according to U.S. Environmental Protection Agency guidelines in effect at

the time for determining virucidal efficacy of disinfectants intended for use on dry inanimate surfaces.

Organic Soil Load: 5% Fetal Bovine Serum Exposure Time: 2 minutes at 68°F

Results: Virucidal against Human Immunodeficiency virus type 1 according to the criteria established by

the U.S. Environmental Protection Agency for registration and labeling of a disinfectant product as a virucide.

**Pathogenic Fungi Efficacy** 

Yeast Organism: Candida albicans [ATCC 14053]

Test Method Used: Modified AOAC Germicidal Spray Method

Organic Soil Load: 5% Horse Serum

Exposure Time: 2 minutes at 72 - 74°F

Incubation: 7 days at 95 - 98.6°F

Results: No growth observed

Yeast Organism: Candida auris AR-BANK#0381 from CDC

Test Method Used: OECD Quantitative Method for Evaluating the Efficacy of Liquid Antimicrobials against Candida auris

on Hard, Non-Porous Surfaces, Wipes and Towelettes

Organic Soil Load: 5% Fetal Bovine Serum Exposure Time: 2 minutes at  $22 \pm 2$  °C Incubation:  $120 \pm 4$  hours at 29 - 31 °C

Results: Met the performance criterion of a minimum reduction in viable cells of 5 Log<sub>10</sub> in accordance with

the U.S. EPA guidance for the Efficacy Evaluation of Products for Claims against *Candida auris*.



### **Toxicity**

### **Acute Inhalation**

Based on the inhalation test results, Super Sani-Cloth Germicidal Disposable Wipe has been classified as Toxicity Category IV for acute inhalation.

### **Acute Oral Toxicity**

Based on the results of this study, Super Sani-Cloth Germicidal Disposable Wipe has been classified as Toxicity Category IV for acute oral toxicity.

#### **Acute Eye Irritation**

Based on the results of this study, Super **Sani-Cloth** Germicidal Disposable Wipe produced eye irritation that indicates the product would be classified as Toxicity Category II for acute eye irritation.

#### **Acute Dermal Toxicity**

Based on the results of this study, Super Sani-Cloth Germicidal Disposable Wipe has been classified as Toxicity Category IV for dermal toxicity.

#### **Acute Dermal Irritation**

Based on the results of primary skin irritation study, Super **Sani-Cloth** Germicidal Disposable Wipe has been classified as Toxicity Category III for dermal effects.

#### **Dermal Sensitization**

Based on the sensitization test results, Super Sani-Cloth Germicidal Disposable Wipe would not be considered a dermal sensitizing agent.

