

The brand you trust, the performance you need.



Ideal for daily use in fast paced environments.

- + Compatible with over 800+ devices and surfaces¹
- + Bactericidal, Tuberculocidal, Virucidal, Fungicidal
- + Intermediate level disinfectant
- + Effective against:
 - Enveloped/non-enveloped viruses, TB, fungi, multi-drug resistant bacteria and bloodborne pathogens^{2,3}
 - *Candida auris*³




DIN 02247354



Super Sani-Cloth®
GERMICIDAL DISPOSABLE WIPE

+PDI®
BE THE DIFFERENCE

+ Color coded for staff identification



The **Dual Access Lid**
Innovation integrated into every PDI canister

+ Enables quick dispensing of a single wipe or multiple wipes



Point of care accessories.



Sani-Bracket®
for wall/pole mounting of a
Sani-Cloth® Canister



Sani-Canister Caddy®
for placement on flat surfaces
(Bilingual)



Wire Softpack Bracket
allows for convenient access



**Educational Wall Charts
and Signage**
(Bilingual)

	REORDER NO.	WIPE SIZE	CASE PACK	CASE WGT	CASE CUBE	PALLET TI/HI
Super Sani-Cloth® Germicidal Disposable Wipe						
Large Canister	P22272	6" X 6.75"	160/can, 12 cans/case	25.63 lbs	1.430 ft	10/3
Extra Large Canister	P33384	7.5" X 15"	75/can, 6 cans/case	15.53 lbs	0.936 ft	10/4



[Learn More](#)

¹Refer to device manufacturer's instructions for use; data on file, ²See Technical Data Bulletin for a complete list of microorganisms, ³On hard, nonporous surfaces

©2025 PDI 140 Woodbridge Avenue, Un CL-1 Market Lane Plaza, Vaughan, ON L4L 4K9
T: 800-263-7067 W: pdihc.ca PDI05259167

