

# Sani-Cloth<sup>®</sup> Bleach

## Ideal for disinfecting high risk areas endemic with Multi-Drug Resistant Organisms (MDRO's), *Clostridium difficile* spores and Norovirus.



Efficacy Claims 50 Disinfection Formulation 1:10 Bleach Dilution

- DIN Registered hospital-grade disinfectant
- Meets CDC and OSHA guidelines
- Bactericidal, Fungicidal, Tuberculocidal, Virucidal
- Compatible with a broad range of surfaces and equipment in healthcare
- Cleans and disinfects in one step

### HealthPRO Contract # CL03843CA

#### DIN 02450410

Bleact

Effective against 50 microorganisms in 4 minutes<sup>1</sup> including<sup>2</sup>

Sani-Cloth

9)(O)(G)

- Norovirus
- Acinetobacter baumannii
- Klebsiella pneumoniae
- Candida albicans
- ESBL Resistant Escherichia coli
- Clostridium difficle
  spores
- MRSA

- Vancomycin Resistant Enterococcus faecalis (VRE)
- Vancomycin Resistant
  Staphyloccocus aureus (VRSA)
- HIV-1<sup>2</sup>
- Hepatitis B virus (HBV)
- Hepatitic C virus (HCV)
- Mycobacterium bovis-BCG (TB)

<sup>1</sup> Contact PDI for additional technical data and refer to the product label for instructions of use.
 <sup>2</sup> Effective against HIV-1 on pre-cleaned environmental surfaces in 60 seconds, when used as directed.





ACCESSIBILITY: Deep well lid provides ample space to store and access next wipe.



**SAFETY:** Icon clearly communicates key safety information on both labels and lids.



**COMPLIANCE:** Prominent contact time for easy identification and improved compliance.

#### DIN 02450410

SANI-CLOTH® BLEACH WIPES	REORDER NO.	PACKING (CASES)	CASE CUBE METERS (FT)	CASE WEIGHT KG (LBS)
Large Canister (15.2 x 26.6 cm)	P24372	12/75's	0.043 (1.54)	10.25 (22.6)
X-Large Canister (19.0 x 38.1 cm)	P24684	12/65's	0.028 (1.01)	9.29 (20.5)
Pail - for Refill (19.0 x 38.1 cm)	P1980P	2/160's	0.03 (0.098)	7.5 (16.5)
<b>Refill</b> (19.0 x 38.1 cm)	P2080P	2/160's	0.15 (0.491)	6.9 (15.2)



## PDI 3-in-1 Universal Sani-Bracket<sup>™</sup>

- Easy to install and place
- Convenient and versatileant
- Promotes usage and compliance

REORDER NO.	PACKING	CASE CUBE	CASE WEIGHT	
	(CASES)	METERS (FT)	KG (LBS)	
P010901	1	0.0048 (0.17)	0.8 (1.75)	





140 Woodbridge Avenue, Unit CL-1, Vaughan, ON L4L 4K9 T: 905-856-4361 T: 1-800-263-7067

©2016 PDI

pdihc.ca PDI06179321