


Infection Prevention Solutions



Our infection prevention solutions target environmental hygiene, hand hygiene and patient care to help get patients back home to their families where they belong.



Environmental Hygiene

PDI Sani-Cloth® Bleach Germicidal Disposable Cloth

- + Bactericidal, fungicidal, tuberculocidal and virucidal
- + Ideal for disinfecting high risk areas endemic with Multi-Drug Resistant Organisms (MDROs), *Clostridium difficile* and Norovirus
- + Effective against the following microorganisms in 4 minutes — Norovirus, *Clostridium difficile* spores, *Acinetobacter baumannii* Multi-Drug Resistant, *Klebsiella pneumoniae* - Carbapenem Resistant, *Klebsiella pneumoniae* - NDM-1 positive, Hepatitis C virus (Human) (HCV) - Bovine Diarrhea Virus, ESBL Resistant *Escherichia coli*, *Mycobacterium bovis*-BCG (TB), and Vancomycin Resistant *Enterococcus faecalis* (VRE)
- + Effective against Human Immunodeficiency Virus Type 1 (HIV-1) on pre-cleaned, hard, nonporous environmental surfaces in 60 seconds, when used as directed
- + Cleans and disinfects



PDI Sani-Cloth® Plus Germicidal Disposable Cloth

- + Bactericidal and virucidal*
- + For use on hard, non-porous surfaces in hospitals and other critical care areas where control of the hazards of cross-contamination is required
- + Effective against the following microorganisms in 5 minutes — *Staphylococcus aureus*, *Salmonella enterica*, *Pseudomonas aeruginosa*, Hepatitis C Virus (HCV), Human Respiratory Syncytial Virus (HRS), Methicillin resistant *Staphylococcus aureus* (MRSA), *Escherichia coli* (E. coli) O157:H7, Influenza A2/Hong Kong, and Vancomycin resistant *Enterococcus faecalis* (VREF)
- + Effective against Human Immunodeficiency Virus Type 1 (HIV-1) on pre-cleaned environmental surfaces in 30 seconds, when used as directed
- + Cleans and disinfects



* Exhibits effective virucidal activity on hard nonporous inanimate surfaces: kills Influenza A2/Hong Kong (representative of the common flu virus) and Herpes Simplex Type II (causative agent of genital herpes) in 5 minutes.

PRODUCT DESCRIPTION	REORDER NO.	PACKING (CASES)	CASE CUBE METERS (FT³)	CASE WGT KG (LBS)
PDI Sani-Cloth® Bleach Germicidal Disposable Cloth				
PDI Sani-Cloth® Bleach Large Canister (15.2 x 26.6cm)	P24372	12/75's	0.0405 (1.43)	11.08 (26.08)
PDI Sani-Cloth® Bleach X-Large Canister (19.0 x 38.1 cm)	P24684	6/65's	0.0265 (0.94)	8.98 (19.79)
PDI Sani-Cloth® Plus Germicidal Disposable Cloth				
PDI Sani-Cloth® Plus Large Canister (15.24 x 17.15 cm)	Q90172	12/160's	0.0405 (1.43)	11.63 (25.63)
PDI Sani-Cloth® Plus X-Large Canister (19.05 x 38.10 cm)	Q85084	6/65's	0.0265 (0.94)	6.71 (14.79)

Super Sani-Cloth® Germicidal Disposable Wipe

- + Bactericidal, tuberculocidal and virucidal*
- + For use in fast-paced environments that require a short contact time and broad coverage of microorganisms
- + Cleans and disinfects
- + Effective against the following microorganisms in 2 minutes — *Staphylococcus aureus*, *Salmonella enterica*, *Pseudomonas aeruginosa*, *Acinetobacter baumannii*, *Klebsiella pneumonia*, *Mycobacterium bovis*-BCG (TB), Hepatitis C virus (Human) (HCV) – Bovine Diarrhea Virus, Methicillin resistant *Staphylococcus aureus* (MRSA), *Escherichia coli* (E. coli) O157:H7, Influenza A2/ Hong Kong (Hong Kong flu), Vancomycin resistant *Enterococcus faecalis* (VRE), Measles Virus and *Candida auris*.



Easy Screen Cleaning Wipe

- + Features a 70% isopropyl alcohol (IPA) formulation
- + Compatible with touchscreen healthcare equipment, tablets and smart phones¹
- + Fast-drying and anti-streaking
- + Cleans dirt, grime, fingerprints, and smudges from touchscreen surfaces¹
- + Meets AORN's "Guideline for Surgical Attire" to clean cell phones, tablets, and other personal communication or hand-held electronic equipment²



* Kills viruses such as Influenza A2 /Hong Kong (Hong Kong flu) and Herpes Simplex Type 2 in 2 minutes

¹ Cleaning the screen of your mobile device with a wipe containing alcohol may reduce the life span of any anti-fingerprint coating. Please check your device manufacturer's cleaning instructions.

² AORN Guideline For Surgical Attire, Section 9. Personal Items, Recommendation 9.2, page 12

PRODUCT DESCRIPTION	REORDER NO.	PACKING (CASES)	CASE CUBE METERS (FT ³)	CASE WGT KG (LBS)
Super Sani-Cloth® Germicidal Disposable Wipes				
Large Canister (6" x 6.75")	P22272	12/160's	0.0405 (1.43)	11.63 (25.63)
X-Large Canister (7.5" x 15")	P33384	6/65's	0.0265 (0.94)	6.30 (13.80)
PDI Easy Screen Cleaning Wipe				
Medium Canister (6" X 9")	P07872	12/70's	0.028 (0.99)	4.90 (10.81)

Skin Antiseptics

PDI Povidone-Iodine Prep Pads and Swabsticks

- + Saturated with a 10% USP Povidone-Iodine prep solution
- + Antiseptic wound cleanser⁴
- + Helps to reduce the risk of infection in minor cuts and scratches⁴



⁴Data on file (C0313)

PRODUCT DESCRIPTION	REORDER NO.	PACKING (CASES)	CASE CUBE METERS (FT ³)	CASE WGT KG (LBS)
PDI POVIDONE-IODINE Prep Pads and Swabsticks				
PDI Povidone-Iodine Prep Pad · Medium	B51200	10/100's	0.0092 (0.325)	2.25 (4.96)
PDI Povidone-Iodine Swabstick 1's	S42050	10/50's	0.3820 (13.49)	3.15 (6.94)
PDI Povidone-Iodine Prep Pad Swabstick 3's	S42325	10/25's	0.4392 (15.51)	3.45 (7.61)

Patient Care

PDI Sani-Hands® Instant Hand Sanitizing Wipes

- + Cleans, sanitizes and moisturizes⁵

Hygea® Multi-Purpose Premium Washcloth

- + Suitable to use on all areas of the body, including face and neck
- + Ideal for incontinence care
- + Gentle cleansing formula contains soothing aloe and vitamin E
- + Hypoallergenic, alcohol-free and not made with natural rubber latex⁶
- + Removes dirt, sweat, soot, and grime from skin



⁵Data on file (C0207)

⁶Data on file (C0113)

PRODUCT DESCRIPTION	REORDER NO.	PACKING (CASES)	CASE CUBE METERS (FT ³)	CASE WGT KG (LBS)
PDI Sani-Hands® Instant Hand Sanitizing Wipes				
PDI Sani-Hands® Instant Hand Sanitizing Wipes 135ct	P12872	12/135's	0.0280 (0.99)	8.3 (18.4)
PDI Sani-Hands® Instant Hand Sanitizing Wipes, Individual Packet	D72400	10/100's	0.0184 (0.65)	4.5 (9.9)
PDI Hygea® Adult Care				
PDI Hygea® Premium Multipurpose Washcloth 60ct	J14143	6/60's	0.0190 (0.67)	5.1 (11.2)

Nice 'N Clean® Baby Wipes

- + Developed for baby's soft and delicate skin⁷
- + Wipe is pre-moistened with a mild, cleansing solution⁷
- + Unscented, alcohol-free*, and does not contain natural rubber latex⁷
- + Features a gentle formula with aloe and vitamin E⁷



*Does not contain ethanol, isopropanol or rubbing alcohol
⁷Data on file (C0287)

PRODUCT DESCRIPTION	REORDER NO.	PACKING (CASES)	CASE CUBE METERS (FT ³)	CASE WGT KG (LBS)
Nice 'N Clean® Baby Care				
PDI Nice 'n Clean® Baby Wipes - Unscented 80ct	M233XT	12/80's	0.0212 (0.75)	5.8 (12.7)
PDI Nice 'n Clean® Baby Wipes - Unscented 40ct	Q70040	12/40's	0.0108 (0.38)	2.9 (6.4)

Specialty Products



PDI Adhesive Tape Remover Pads



PDI See Clear Lens Cleaning Wipes



PDI Castile Soap Towelettes

PRODUCT DESCRIPTION	REORDER NO.	PACKING (CASES)	CASE CUBE METERS (FT ³)	CASE WGT KG (LBS)
PDI Adhesive Tape Remover Pads				
Individual Packet (1.2 x 2.6")	B16400	10/100's	0.0071 (0.25)	1.6 (3.6)
PDI See Clear Lens Cleaning Wipes				
Individual Packet	D25431	12/120's	0.0252 (0.89)	4.4 (9.6)
PDI Castile Soap Towelettes				
Individual Packet	D41900	10/100's	0.0184 (0.65)	5.4 (12.0)

Point-of-Care Accessories

- + Bracket can be used for Sani-Hands® or Sani-Cloth® Plus products
- + Easy to install and place
- + Convenient and versatile
- + Promotes usage and compliance



PDI helps reduce preventable infections, control healthcare costs, and ultimately help save lives by delivering a broad range of evidence-based, market leading environmental, hand hygiene and patient care solutions. Our expertise in making a difference in the community and healthcare environments extends worldwide.

For more information, or to place an order,
contact your local PDI sales representative.

PDI®, Sani-Cloth®, Sani-Hands®, Hygea®, Be The Difference®, Sani-Bracket™ and See Clear, and all associated logos and product images, are trademarks or properties of Professional Disposables International, Inc. Nice 'N CLEAN® is a registered trademark of Nice-Pak Products, Inc.

©2025 PDI 0925 UPDATE PDI05158054
140 Woodbridge Avenue, UN CL-1 Market Lane Plaza Vaughan (Ontario) L4L 4K9
Homologation ISO 9001:2008 + homologation ISO 13485:2003
T. 905 856-4361 + t 800 263-7067 (au Canada seulement)
www.pdihc.ca

